



**\*22852\***

**22852**

PATENT TRADEMARK OFFICE

**RECEIVED**  
APR 08 2003  
TECH CENTER 160012900

PATENT

Attorney Docket No. 03806.0517

**IN THE UNITED STATES PATENT AND TRADEMARK OFFICE**

In re Application of:

Francis BLANCHE *et al.*

Application No.: 09/970,663

Filed: October 5, 2001

For: COMPOSITION FOR THE  
PRESERVATION OF  
INFECTIOUS RECOMBINANT  
ADENOVIRUSES

Group Art Unit: 1635

Examiner: Brian Whiteman

Assistant Commissioner for Patents  
Washington, DC 20231

Sir:

**DECLARATION UNDER 37 C.F.R. § 1.131**

I, Francis Blanche, state that I am one of the named applicants of the above-identified application and that I am a co-inventor of the subject matter described and claimed therein. Prior to November 16, 1998, we, the co-inventors had completed in France the invention as described and claimed in the above-identified application as evidenced by the following:

1. Exhibit A: Laboratory Notebook Pages 51-55 and 176 (A1-A6) of Francis Blanche, showing, a composition comprising adenoviral particles and a glycerol buffer solution at pH 8.4, wherein the buffer solution does not contain added divalent metal cations or alkali metal cations. See pages 52-53 (A2-A3), formulation #2, for example, comprises Tris/HCl and 10% glycerol at pH 8.4 (hereinafter referred to as "formulation #2".) The addition of adjuvants, such as sucrose or Tween20 is shown, for example, at page 176, formulations C and D. Formulation #2 is shown to be

useful for preserving adenoviruses. See page 55 (A5), stable viral titer after 15 days of storage in formulation #2. Some compositions were tested for stability after  $-20^{\circ}\text{C}$  or  $4^{\circ}\text{C}$  storage, indicating that the  $-20^{\circ}\text{C}$  frozen viral compositions were thawed to test viability. See page 176 (A6), last three lines from the bottom.

2. The present specification at page 17, first formulation in the Table, shows a formulation identical to formulation #2 of Exhibit A;
3. Example 3 of the present specification, at pages 18-19, shows that a formulation identical to formulation #2 of Exhibit A has a stable viral titer after 15 days of storage, similar to the 15-day storage stability of formulation #2 shown on page 55 (A5) of Exhibit A.

While the dates have been redacted, the undersigned testifies that all experiments described herein were conducted before November 16, 1998.

I declare further that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further, that the statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patents issuing thereon.

Dated: March 28, 2003

By: \_\_\_\_\_

  
Francis Blanche

## EXHIBIT "A" (6 pages)

051

ESSAI N° \_\_\_\_\_

GEL 02051

ESSAIS FORMULATIONS  
STABILITE

BUT : Observer la stabilité ou la précipitation éventuelle du virus Y28 dans différentes formulations.

MATERIEL VIRAL ETUDIE :

Solution virale Y28 produite en Cell Cube à l'échelle 8 Mer par l'équipe JF Chaubard et purifiée par chromatographie échangeuse d'anions, conservée dans le Tris 20mM pH8, MgCl<sub>2</sub> 1mM, NaCl 500mM et glycerol 10%.

Le virus purifié titre  $3,94 \cdot 10^{11}$  pv/ml.

PREPARATION DES DIFFERENTS TAMPONS ETUDIES :1. Solutions mères :

	SOLUTIONS MERES :	PREPARATIONS :
A	Tris / HCl pH 8,4 à 500mM	10,07g Tris base + 6,60g Tris/HCl dans 250ml eau PPI (Tris base ref: T8524 et Tris HCL ref: T7149)
B	Sucrose à 50g/100ml	250g de sucrose dans 500ml d'eau PPI
C	NaCl 5M	Sigma - Aldrich ref: S150
D	MgCl <sub>2</sub> 1M	Sigma - Aldrich ref: M1028
E	Glycerol	Sigma - Aldrich ref: G5516
F	D-Mannitol	Sigma - Aldrich ref: M9647
G	Tween 20	Sigma - Aldrich ref: P8074
H	Tampon borate pH 7.4 100mM	Acide borique 100mM + NaOH 0,1N
I	Tampon phosphate pH 7.4 10mM	130mg KH <sub>2</sub> PO <sub>4</sub> + 705mg K <sub>2</sub> HPO <sub>4</sub> dans 500ml eau PPI

A-1

052

ESSAI N° \_\_\_\_\_

B1

2. Formulations :

SOLUTIONS MÈRES :										
	A	B	C	D	E	F	G	H	I	Eau PPI
ESSAIS :										gsp 500ml
1	20ml									gsp 500ml
2	20ml				+ 50ml					gsp 500ml
3	20ml	50ml		0,5ml						gsp 500ml
4	20ml	50ml								gsp 500ml
5	20ml	50ml		0,5ml		25g				gsp 500ml
6	20ml	50ml	15ml	0,5ml		25g				gsp 500ml
7	20ml	50ml					0,5ml			gsp 500ml
8	20ml	50ml		0,5ml			0,5ml			gsp 500ml
9		50ml		0,5ml				50ml		gsp 500ml
10					+ 50ml				500ml	
11	Solution virale obtenue au 2ème rinçage lors de la dialyse finale. Titre = $2,88.10^{11}$ pv/ml dans DPBS/ NaCl 150mM / glycérol 10%									

3 Résumé des formulations étudiées :

Voir tableau ci-après.

053

ESSAI N°

ESSAIS	Tch 20mm pH 8.4	NaCl 150mM	MeCl 1mM	Sucrose 5%	Tween20 0.1%	Mannitol 5%	Glycerol 10%	DPBS	Tampon libre borate pH 7.4	Tampon libre pH 7.4 phosphate
1	+									
2	+						+			
3	+			+						
4	+					+				
5	+		+	+						
6	+			+						
7	+			+					+	
8	+			+						+
9										
10										
11		+								

V28 CELL CUBE BMEK

ESSAIS DE FORMULATIONS

A-3

054

SAI N° \_\_\_\_\_

MATÉRIELS UTILISÉS :

- 10 PD 10 pour la diafiltration équilibrée avec 5 x 5 ml de tampon étudié.
- Ultrafree 15 ml avec membrane Biomax 100 Kd (Millipore) (2x pour chaque essai).
- Centrifugeuse réglée à 1500  $\frac{r}{mn}$ .

MISE EN ŒUVRE :

<u>OPÉRATIONS :</u>	<u>POUR CHAQUE ESSAI :</u>
<u>DIAFILTRATION :</u>	10 PD10 x 2,5 ml de solution virale Y28 à $3,94 \cdot 10^{11}$ pv/ml. Elution par 10 PD10 x 2,5 ml du tampon étudié.
<u>CONCENTRATION :</u>	2 Ultrafree 15 ml 100Kd remplie à 15 ml puis rechargés avec 2,5 ml de solution virale diafiltrée. Soit 17,5 ml concentrés à 500 $\mu$ l (x2). (soit une concentration à $\approx 1 \cdot 10^{12}$ pv/ml.)
<u>RÉCUPÉRATION ET FILTRATION 0,2 <math>\mu</math>m :</u>	Récupération en pool des 2 Ultrafree pour chaque essai. Filtration au filtre Millex 0,2 $\mu$ non stériles. Stockage dans tubes en verre stériles.
<u>ALIQUOTAGE : (t=0)</u>	→ 100 $\mu$ l dans tube Eppendorf congelé à -26°C par essai. → 20 $\mu$ l + 980 $\mu$ l tampon chip anal. pour dosage. → env. 900 $\mu$ l conservés à +4°C pour étude de stabilité. → env. 100 $\mu$ l de la solution virale Y28 sortie chromatographique initiale est congelé à -26°C.
<u>TEMOINS PBS/tycerno 10% :</u>	directement concentré à $1 \cdot 10^{12}$ pv/ml, récupéré et aliquoté comme les autres essais.

A-4

055

ESSAI N°

## DOSAGES CLAP ANALYTIQUE

ISSAISTAMPONS	TITRE $\mu\text{V/ml } J=0$	APPARENCE DE L'ÉCHANTILLON	TITRE $\mu\text{V/ml } J=15$	OBSERVATIONS	TITRE $\mu\text{V/ml } J=20$ (c)	TITRE $\mu\text{V/ml } J=22$ (c)
Tampou 1 Tch 20mm	4.97.10 <sup>11</sup>	normale à J=15	non filtré : 1.0.10 <sup>11</sup> filtré 0.2µm : 9.1.10 <sup>11</sup>	le retour pic adeno en face nive plateau : 2000 synthétique : 25 et 1.5	non dosé normal	—
Tampou 2 Tch 20mm	7.71.10 <sup>11</sup>	normale à J=15	non filtré : 8.12.10 <sup>11</sup> filtré 0.2µm : 7.96.10 <sup>11</sup>	pic synthétique	non filtré : 7.88.10 <sup>11</sup> pic synthétique normal	non filtré : 9.27.10 <sup>11</sup> pic synthétique normal
Tampou 3 Tch 20mm	6.29.10 <sup>11</sup>	opacification à J=12 <sup>1</sup> mais non précipité à J=15	non filtré : 2.31.10 <sup>11</sup> filtré 0.2µm : 3.09.10 <sup>11</sup>	normal du pic synthétique le retour plateau : 2000 synthétique : 0.93 et 0.86	non dosé normal mais non J	—
Tampou 4 Tch 20mm	6.31.10 <sup>11</sup>	normale à J=15	non filtré : 5.83.10 <sup>11</sup> filtré 0.2µm : 5.7.10 <sup>11</sup>	pic synthétique	non filtré : 1.87.10 <sup>11</sup> nive plateau : 14000 synthétique : 28 et 1.12	non filtré : 1.09.10 <sup>11</sup> nive plateau : 4000 synthétique : 0.87 et 0.68
Tampou 5 Tch 20mm	5.84.10 <sup>11</sup>	normale à J=15	non filtré : 1.85.10 <sup>11</sup> filtré 0.2µm : 1.47.10 <sup>11</sup>	le retour pic adeno en face nive plateau : 17000 synthétique : 1.08 et 1.12	non dosé normal	—
Tampou 6 Tch 20mm	6.48.10 <sup>11</sup>	précipité à J=7	non filtré : non dosé filtré 0.2µm : —	—	—	—
Tampou 7 Tch 20mm	6.22.10 <sup>11</sup>	normale à J=15	non filtré : 9.53.10 <sup>11</sup> filtré 0.2µm : 9.21.10 <sup>11</sup>	semaine du pic adeno nive plateau : 9000 synthétique : 0.95 et 0.83	non dosé normal	—
Tampou 8 Tch 20mm	7.17.10 <sup>11</sup>	opacification à J=7 précipité le lendemain	non filtré : non dosé filtré 0.2µm : —	—	—	—
Tampou 9 Tch 20mm	non dosé	virus retenus sur le filtre solution insoluble des le changement de l'ampou	non filtré : non dosé filtré 0.2µm : —	—	—	—
Tampou 10 Tch 20mm	7.20.10 <sup>11</sup>	opacification à J=2 précipité le lendemain	non filtré : non dosé filtré 0.2µm : —	—	—	—
Tampou 11 Tch 20mm	5.37.10 <sup>11</sup>	précipité à J=1 jour	non filtré : non dosé filtré 0.2µm : —	—	—	—

note: pour le pic adeno en face - même plateau : 3000 / synthétique : 1.1 et 1.16  
(\*) calcul des filtres avec le nouveau étalon J=1  
Essais Tpe et Tpe réalisés à J=22 pour test blancité par M. Laisson

A-5

PAGE 5 of 6

**SUJET :**

## Adenovirus

Echantillon de départ: 400ml fraction F3 (+10% glycérol) du DEMOBATCH 3 (CC16M-Ad5CMV/P53/293), dose  
à  $3,6 \cdot 10^{11}$  pv/ml soit  $1,44 \cdot 10^{11}$  pv pour 400ml.

Tampons étudiés (filres 0,22  $\mu$ m):

- Tampon A : Tris 20mM-pH8,4+10% glycérol
- Tampon B : Tris 20mM-pH8,4+5% sucrose
- Tampon C : Tris 20mM-pH8,4+10% glycérol+5% sucrose
- Tampon D : Tris 20mM-pH8,4+5% glycérol+10% sucrose
- Tampon E : Tris 20mM-pH8,4+10% glycérol+1mM MgCl<sub>2</sub>
- Tampon F : Tris 20mM-pH8,4+10% glycérol+150mM NaCl+1mM MgCl<sub>2</sub>
- Tampon G : Tris 20mM-pH8,4+5% glycérol
- Tampon H : Tris 20mM-pH8,4+10% sucrose

- Tampon I : Acétate d'ammonium 20mM-pH8+10% glycérol  
-Tampon J : Acétate d'ammonium 20mM-pH8+5% sucrose

Mise en place des essais : dans labo L3 de recherches/BI Monod

- Misc en place des échantillons
- 1<sup>ère</sup> étape : concentration de l'échantillon en utilisant 16 Ultrafree 15ml/30Kd membrane biomax (UFV2BTK40 Millipore), centrifugation à 1500n/min. Premier passage, on amène le volume à 5ml, (il faut environ 30mn pour le passage de 15 ml)  
on recharge une deuxième fois les Ultrafree avec 10ml (on tourne à 1760n/min-500G)  
et on amène le volume total final à 105ml.  
on conserve 5ml pour électrophorèse 2D et on effectue un dosage HPLC (d1/10)  
on trouve  $1,21.10^{12}$  pv/ml soit  $1,27.10^{14}$  pv pour 105ml.
  - 2<sup>ème</sup> étape : changement de tampon sur PD10 Pharmacia (4 PD10 par tampon, soit 4 fois 2,5ml du concentrat ou  $1,21.10^{12}$  pv/tampon), on récupère 14ml.
  - 3<sup>ème</sup> étape : on concentre les éluats PD10 sur un Ultrafree 15ml/30Kd (même réf. que étape 1), on amène le volume à <1ml.  
on récupère le concentrat et on volume à 1ml avec le filtrat
  - 4<sup>ème</sup> étape : on fait subir à chaque échantillon une filtration stérilisante sur un filtre Millipore (Sterile Millex-GV 0,22µm) membrane PVDF, récupération dans un tube stérile.
  - 5<sup>ème</sup> étape : sur chaque échantillon de 1ml après filtration → dosage HPLC (d1/50)  
pour les échantillons Tpa à E, aliquotter 14 tubes de 50µl dans tubes stériles,  
pour les échantillons Tpf à J, il y a 15 aliquotes de 50µl.  
les titres se situent entre  $9,8.10^{12}$  et  $1,08.10^{13}$  pv/ml (voir cahier DOS-01 page 42)
  - 6<sup>ème</sup> étape : les aliquotes de 50µl sont mis ce jour en stabilité à -20°C.  
les reliquats soit ~250 à 300µl sont conservés à 4°C.

Il est prévu un dosage pfu (labo D.Faucher) de chaque échantillon → 1 tube de 50µl à -20°C

Nota: Tit. de reactiile a 4°C. pentru cuprinsul noilor.

A-6



**ENGLISH-LANGUAGE TRANSLATION OF  
EXHIBIT "A" (6 pages)**

TRIAL NO. \_\_\_\_\_

CEL 02051

**FORMULATION TRIALS:  
STABILITY.**

**OBJECTIVE:** Observe the stability, or possible precipitation, of the Y28 virus in different formulations.

**VIRAL MATERIAL STUDIED:**

Y28 solution produced in a cell cube on an 8 mer scale by the J.F. Chaubard team, purified by ion exchange chromatography, and preserved in 20mM pH8 TRIS, 1mM MgCl<sub>2</sub>, 500mM NaCl, and 10% glycerol. The purified virus titrates  $3.94 \cdot 10^{11}$  pv/ml.

**PREPARATION OF THE DIFFERENT BUFFER SOLUTIONS USED:**

**1. Stock solutions:**

	STOCK SOLUTIONS:	PREPARATIONS:
A	Tris / HCl pH 8.4 at 500mM	10.07g Tris base + 6.60g Tris/Hcl in 250ml water for injection (Tris base ref: T8524 and Tris HCL ref:T7149)
B	Sucrose at 50g/100ml	250g sucrose in 500ml of water for injection.
C	NaCl 5M	Sigma - Aldrich ref. S150
D	MgCl <sub>2</sub> 1M	Sigma - Aldrich ref. M1028
E	Glycerol	Sigma - Aldrich ref. G5516
F	D-Mannitol	Sigma - Aldrich ref. M9647
G	Tween 20	Sigma - Aldrich ref. P8074
H	100mM borate buffer solution pH 7.4	100mM boric acid + NaOH 0.21N
I	10mM phosphate buffer solution pH 7.4	130mg KH <sub>2</sub> PO <sub>4</sub> + 705mg K <sub>2</sub> HPO <sub>4</sub> in 500ml water for injection.

TRIAL NO. \_\_\_\_\_

2. Formulations:

<u>STOCK SOLUTIONS :</u>										
	A	B	C	D	E	F	G	H	I	Water for injection
TRIAL:										
1	20ml									QS 500ml
2	20ml				+ 50ml					QS 500ml
3	20ml	50ml		0,5ml						QS 500ml
4	20ml	50ml								QS 500ml
5	20ml	50ml		0.5ml		25g				QS 500ml
6	20ml	50ml	15ml	0.5ml		25g				QS 500ml
7	20ml	50ml					0.5ml			QS 500ml
8	20ml	50ml		0.5ml			0.5ml			QS 500ml
9		50ml		0.5ml				50ml		QS 500ml
10					+ 50ml				500ml	
11	<u>Viral solution obtained in the second rinsing during the final diafiltration.</u> <u>Titer = <math>2.8810^{11}</math> pv/ml in DPBS/1500mM NaCl/glycerol 10%.</u>									

3. Summary of the formulations studied:

See the following tables.

TRIAL NO. \_\_\_\_\_

**Y28 CELL CUBE 8MER****FORMULATION TRIALS**

<u>TRIAL</u>	<u>Tris 20mM pH 8(?)</u>	<u>NaCl 150mM</u>	<u>MgCl<sub>2</sub> 1mM</u>	<u>Sucrose 5%</u>	<u>Tween20 0.1%</u>	<u>Mannitol 5%</u>	<u>Glycerol 10%</u>	<u>DPBS</u>	<u>10 mM borate pH7.4 buffer solution</u>	<u>10mM pH7.4 phosphate buffer solution</u>
1	+									
2	+						+			
3	+		+	+						
4	+			+						
5	+		+	+		+				
6	+	+	+	+		+				
7	+			+	+					
8	+		+	+	+					
9			+	+					+	
10							+			+
11		+					+	+		

TRIAL NO. \_\_\_\_\_

MATERIALS USED:

- 10 PD 10 for diafiltration balanced with 5 x 5ml of the buffer solution studied.
- 15ml Ultrafree with 100 Kd Biomax (Millipore) membrane (2x for each trial).
- Centrifuge set at 1500 rev/min.

IMPLEMENTATION:

<u>OPERATIONS:</u>	<u>FOR EACH TRIAL:</u>
<u>DIAFILTRATION:</u>	10 PD10 x 2.5ml of Y28 viral solution at $3.94 \cdot 10^{11}$ pv/ml. Elution by 10PD10 x 3.5ml of the buffer solution studied.
<u>CONCENTRATION:</u>	15ml 100Kd 2 Ultrafree filled to 15ml and then refilled with 2.5 diafiltrated viral solution. 17.5ml concentrated at 500 $\mu$ l (x2). (or a concentration at $\approx 1 \cdot 10^{13}$ pv/ml.)
<u>RECOVERY AND FILTRATION 0.2<math>\mu</math>m:</u>	Recovery and pooling of the 2 Ultrafree for each trial. Filtration using unsterilized 0.2 $\mu$ Millex filters. Storage in sterilized glass tubes.
<u>ALIUQUOTING : (t=0)</u>	→ 100 $\mu$ l in Ependorff tube frozen at -26°C. → * → 20 $\mu$ l + 980 $\mu$ l anal. HPCL buffer solution for dosing. → About 900 $\mu$ l stored at +4°C to study stability. → About 100 $\mu$ l of the initial chromate emerging Y28 viral is frozen at -26°C.
<u>10% glycerol/PBS Samples:</u>	Frozen directly at $1 \cdot 10^{13}$ pv/ml, recovered and aliquoted in the same way as the other trials.

TRIAL NO. \_\_\_\_\_

ANALYTICAL HPLC MEASUREMENTS:

<u>TRIALS/BUFFERS</u>	<u>TITER pV/ml J=0</u>	<u>SAMPLE APPEARS</u>	<u>TITER pV/ml day=15</u>	<u>OBSERVATIONS HPLC of the dosage, day=15</u>	<u>TITER pV/ml day=20(*) HPLC Observations Sample appears</u>	<u>TITER pV/ml day=22(*) HPLC Observations Sample appears</u>
<u>Buffer 1</u> Tris 20mM	4.97. 10 <sup>12</sup>	normal at day=15	unfiltered: 1.0. 10 <sup>12</sup> filtered 0.2µm: 9.1. 10 <sup>11</sup>	The adeno return peak trails plate number: 12,000 asymmetries: 1.25 and 1.50	undosed normal	—
<u>Buffer 2</u> Tris*glycerol	7.71. 10 <sup>12</sup>	normal at day=15	unfiltered: 8.12. 10 <sup>12</sup> filtered 0.2µm: 7.96. 10 <sup>12</sup>	symmetrical peak	unfiltered: 7.88. 10 <sup>13</sup> normal symmetrical peak	unfiltered: 9.27. 10 <sup>12</sup> normal symmetrical peak
<u>Buffer 3</u> Tris+MgCl <sub>3</sub> + sucrose	6.29. 10 <sup>12</sup>	opacification at day=12 <sup>3</sup> but not precipitated at day=15	unfiltered: 2.33. 10 <sup>11</sup> filtered 0.2µm: 2.09. 10 <sup>11</sup>	asymmetrical rise of the peak plate number: 32,000 asymmetries: 0.93 at 0.86	undosed clouding but not i	—
<u>Buffer 4</u> Tris+sucrose	6.31. 10 <sup>12</sup>	normal at day=15	unfiltered: 5.83. 10 <sup>12</sup> filtered 0.2µm: 5.7. 10 <sup>12</sup>	symmetrical peak	unfiltered: 1.87. 10 <sup>12</sup> plate number: 14,000 asymmetries: 1.28 and 1.42 (normal)	unfiltered: 1.09. 10 <sup>12</sup> plate number: 4,000 asymmetries: 0.87 and 0.68 (normal)
<u>Buffer 5</u> Tris+MgCl <sub>2</sub> +sucrose+mammitol	5.84. 10 <sup>12</sup>	normal at day=15	unfiltered: 1.85. 10 <sup>12</sup> filtered 0.2µm: 1.47. 10 <sup>12</sup>	The adeno return peak trails plate number: 17,000 asymmetries: 1.08 and 1.12	undosed normal	—
<u>Buffer 6</u> Tris+NaCl+MgCl <sub>2</sub> +sucrose+mammitol	6.48. 10 <sup>12</sup>	precipitated at day=7	unfiltered: undosed filtered 0.2µm:	—	—	—
<u>Buffer 7</u> Tris+sucrose+Tween	6.22. 10 <sup>12</sup>	normal at day=15	unfiltered: 9.53. 10 <sup>11</sup> filtered 0.2µm: 9.31. 10 <sup>11</sup>	rounded peak top plate number: 9,000 asymmetries: 0.95 and 0.83	undosed normal	—
<u>Buffer 8</u> Tris+MgCl <sub>2</sub> +sucrose+Tween	7.17. 10 <sup>12</sup>	opacification at day=7 precipitated the next day.	unfiltered: undosed filtered 0.2µm:	—	—	—
<u>Buffer 9</u> Borate+MgCl <sub>2</sub> +sucrose	Undetected	virus held on the filter solution clouds once the buffer solution is changed.	unfiltered: undosed filtered 0.2µm:	—	—	—
<u>Buffer 10</u> Phosphate + glycerol	7.20. 10 <sup>12</sup>	opacification at day=2 precipitated the next day	unfiltered: undosed filtered 0.2µm:	—	—	—
<u>Buffer 11</u> DPBS+NaCl+glycerol	5.37. 10 <sup>12</sup>	precipitated at < 1 day	unfiltered: undosed filtered 0.2µm:	—	—	—

Note: for the adeno return peak measurement standard → plate number 32,000/asymmetries: 1.1 and 1.16.

(\*) computation of titers with the new measurement standard: 141

TRIAL NO. \_\_\_\_\_

SUBJECT: ADENOVIRUS**CONDUCTING ADENOVIRUS STABILITY TRIALS IN DIFFERENT FORMULATIONS**

Starting sample: 400ml fraction F3 (+10% glycerol) of DEMOBATCH 3 (CC16M-Ad5/CMV/P53/293) dosed at  $3.6 \cdot 10^{11}$  pv/ml or  $1.44 \cdot 10^{14}$  pv per 400ml.

Buffer solutions studied (0.22µm filtered):

- Buffer solution A: Tris 20mM-pH8.4+10% glycerol
- Buffer solution B: Tris 20mM-pH8.4+5% sucrose
- Buffer solution C: Tris 20mM-pH8.4+10% glycerol+5% sucrose
- Buffer solution D: Tris 20mM-pH8.4+5% glycerol+10% sucrose
- Buffer solution E: Tris 20mM-pH8.4+10% glycerol+1mM  $MgCl_2$
- Buffer solution F: Tris 20mM-pH8.4+10% glycerol+150mM NaCl+1mM  $MgCl_2$
- Buffer solution G: Tris 20mM-pH8.4+5% glycerol
- Buffer solution H: Tris 20mM-pH8.4+10% sucrose

- Buffer solution 1: ammonium acetate 20mM-pH8+10% glycerol
- Buffer solution 1: ammonium acetate 20mM-pH8+5% sucrose

Carrying Out the Trials: At Research Lab L3/Bt Monod

1<sup>st</sup> Step: Concentrating the sample by using 15ml/30Kd 16 Ultrafree biomax membrane (UFV2BTK40 Millipore), centrifuged at 1500rev/min. First run, volume brought to 5ml (5ml run requires @30 mins). The Ultrafree is filled a second time with 10ml (turning occurs at 1760 rev/min.-500G). The final total volume is brought to 105ml. 5ml is stored for 2D electrophoresis and HPLC (dl/10) measurement occurs. One then finds  $1.21 \cdot 10^{12}$  pv/ml, or  $1.27 \cdot 10^{14}$  pv per 105ml.

2<sup>nd</sup> Step: Changing over the sample to PD10 Pharmacia (4 PD10 by buffer solution, i.e., 4 x 2.5ml of the concentrate or  $1.21 \cdot 10^{13}$  pv/buffer solution), 14ml are recovered.

3<sup>rd</sup> Step: The PD10 eluates are concentrated on a 15ml/30Kd Ultrafree (same ref. as Step 1) and the volume is brought to <1ml. The concentrate is recovered and the volume is increased to 1ml with filtrate.

4<sup>th</sup> Step: Each sample undergoes a sterilizing filtration on a Millipore film (Sterile Millex-GV 0.22µm) membrane (PVDF). Collected in a sterile tube.

5<sup>th</sup> Step: On each 1ml sample after filtration →HPLC (dl/50). For samples TpA to E, aliquot 14 tubes of 50µl in sterile tubes. For samples TpF to J, there are 15 aliquots of 50µl. The titers are located between  $9.8 \cdot 10^{12}$  and  $1.08 \cdot 10^{13}$  pv/ml (see Manual DOS-01 page 42).

6<sup>th</sup> Step: The 50µm aliquots are used while stable at -20°C. The carry-over, i.e., 250 to 300µl, is stored at 4°C.

A PFU (D. Faucher Lab) measurement of each sample is provided →1 tube of 50µl at -20°C.